

Lab 6

Complexometric Determination

Partner: Linh Nguyen

45

10/22/2019

Prelab Assignment

Concentration of EDTA,  $[EDTA] = \frac{3.7145g}{0.1000L} \left( \frac{1mol}{372.24g} \right) = 9.979 \times 10^{-3} M$

At pH 10, EDTA reacts with both  $Mg^{2+}$  &  $Ca^{2+}$

At pH 12, EDTA only reacts with  $Ca^{2+}$

pH 10

Titrated volume EDTA = 11.23 mL = 0.01123 L

mol EDTA =  $0.01123L \left( \frac{9.979 \times 10^{-3} mol}{1L} \right) = 1.121 \times 10^{-4} mol$

Since the reaction is 1:1 ratio (between EDTA &  $Mg^{2+}$  &  $Ca^{2+}$ ),

mol of  $Mg^{2+}$  + mol of  $Ca^{2+}$  =  $1.121 \times 10^{-4} mol$

pH 12

Titrated volume = 2.11 mL = 0.00211 L

mol EDTA =  $0.00211L \left( \frac{9.979 \times 10^{-3} mol}{1L} \right) = 2.11 \times 10^{-5} mol$

Since the reaction is 1:1 ratio (between EDTA &  $Ca^{2+}$ ),

mol of  $Ca^{2+}$  =  $2.11 \times 10^{-5} mol$

mol of  $Mg^{2+}$  =  $1.121 \times 10^{-4} mol - 2.11 \times 10^{-5} mol = 9.10 \times 10^{-5} mol$

$[Mg^{2+}] = [MgSO_4] = \frac{9.10 \times 10^{-5} mol}{0.100L} = 9.10 \times 10^{-4} M \cdot \left( \frac{120.367g}{1mol} \right) \left( \frac{1000mg}{1g} \right) = 110mg/L$

$[Ca^{2+}] = [CaCl_2] = \frac{2.11 \times 10^{-5} mol}{0.100L} = 2.11 \times 10^{-4} M \cdot \left( \frac{110.984g}{1mol} \right) \left( \frac{1000mg}{1g} \right) = 23.4mg/L$





46

Lab 6/2019

Partner: Linh Nguyen

## Compleximetric Determination

10/22/2019

Mass of EDTA (g)  $\rightarrow$  0.8653Mass of empty beaker (g)  $\rightarrow$  119.8945Mass of beaker + EDTA (g)  $\rightarrow$  120.7598Molarity of EDTA solution (M)  $\rightarrow$  0.009298



Complexometric Determination

10/22/2019

Objective

The purpose of this experiment is to determine the concentration of calcium ions and magnesium ions in a water sample.

Introduction

The concentration of magnesium and calcium ions in the water sample can be determined by titrimetrically mixing a known concentration of EDTA into a ~~set~~ the water sample containing an indicator. This determination is done in two parts, where

Eriochrome Black T is used as the indicator in the first part and hydroxynaphthol blue is used as the indicator in the second part.

The first part will determine the sum of <sup>concentration of</sup> calcium and magnesium ions in the water sample at pH 10, where the titrant turned blue from red at end point. The second part will determine only the concentration of calcium ions in the water sample, as the magnesium ions have reacted with  $\text{OH}^-$  to form solid, at pH 12, where the titrant ~~turned~~ retained the blue color for a ~~are~~ long period. From there, the magnesium ions can be calculated from two measured/calculated values.

Procedure

1. Collect a 50mL water sample in a clean, dry beaker.
2. Weigh about 0.9g of sodium EDTA solid to the nearest 0.1mg. & place the weighed EDTA in a 250mL beaker.
3. Dissolve the EDTA in about 100mL of distilled water, which will occur slowly but can be speeded up by warming the solution gently.



48

Lab 6

Partner: Linh Nguyen

Complexometric Determination

10/22/2019

Piscicola

Eriochrome Black T as indicator

Sample	Reference	1	2	3
Initial	/			
buret reading (mL)		0.31	0.69	13.35
Final buret reading (mL)		12.97	13.35	26.02
Titrated volume (mL)		12.66	12.66	12.67

Eriochrome Blank T indicator			
Sample number	1	2	3
Volume of sample water (mL)	75.0	75.0	75.0
Initial buret reading (mL)	0.31	0.69	13.35
Final buret reading (mL)	12.97	13.35	26.02
Titrated volume (mL)	12.66	12.66	12.67
Mol EDTA titrated (mol)	1.177E-04	1.177E-04	1.178E-04
Total mol of $Mg^{2+}$ and $Ca^{2+}$ (mol)	1.177E-04	1.177E-04	1.178E-04





## Compleximetric Determination

10/22/2019

4. After the EDTA dissolves, transfer ~~the~~ the solution quantitatively to a 250mL volumetric flask. Fill the flask to the mark with water, insert the stopper & shake to mix the solution thoroughly.
5. Make a reference sample containing 75.0mL of DI water, the pH 10 buffer (2mL), & the indicator (2-4 drops of Eriochrome Black T). Use this sample to find out the end point of titration.
6. Use a pipet (3 x 25.0mL) to deliver 75.0mL water sample into each of 3, 250mL Erlenmeyer flasks.
7. Add 2mL of pH 10 buffer & 2-4 drops of Eriochrome Black T indicator solution to each flask (solutions should be <sup>red</sup>).
8. Fill a 50mL buret with the EDTA solution & record the initial volume to the nearest 0.01mL. Titrate each water sample to the end point where the titrant solution changes abruptly from red to blue. Record each end-point volume to the nearest 0.01mL.
9. Make a reference sample containing 75.0mL of DI water, the NaOH solution (10mL of 0.1M), & the indicator (0.1g of calcon). Use this sample to find out the end point of titration.
10. Again pipet 75.0mL of water into each of 3 250mL Erlenmeyer flasks. Add 10.0mL of 0.1M NaOH into each flask using a graduated cylinder, and swirl to mix well. A white precipitate of  $Mg(OH)_2$  may be visible at this point.



S.O

Lab 6

Compleximetric determination

10/22/2019

Partner = Linh Nguyen

Calcon as indicator

Sample	Reference	1	2	3
Initial buret reading (mL)	/	0.68	11.4	22.19
Final buret reading (mL)		11.4	22.19	32.95
Titrated volume (mL)		10.72	10.79	10.76

Calcon Indicator			
Sample number	1	2	3
Volume of sample water (mL)	75.0	75.0	75.0
Initial buret reading (mL)	0.68	11.4	22.19
Final buret reading (mL)	11.4	22.19	32.95
Titrated volume (mL)	10.72	10.79	10.76
Mol EDTA titrated (mol)	9.968E-05	1.003E-04	1.000E-04
Mol of $\text{Ca}^{2+}$ (mol)	9.968E-05	1.003E-04	1.000E-04
Mol of $\text{Mg}^{2+}$ (mol)	1.807E-05	1.742E-05	1.770E-05
Mass of $\text{CaCl}_2$ (mg)	1.106E+01	1.113E+01	1.110E+01
Mass of $\text{MgSO}_4$ (mg)	2.1750	2.0966	2.1302
Concentration $\text{CaCl}_2$ in water (mg/L)	147.5	148.5	148.0
Average concentration $\text{CaCl}_2$ (mg/L)	148.0		
Standard deviation of $\text{CaCl}_2$ (mg/L)	0.5		
RSD of $\text{CaCl}_2$ (%)	0.33		
Concentration $\text{MgSO}_4$ in water (mg/L)	29.0	28.0	28.4
Average concentration $\text{MgSO}_4$ (mg/L)	28.5		
Standard deviation of $\text{MgSO}_4$ (mg/L)	0.5		
RSD of $\text{MgSO}_4$ (%)	1.8		





## Compleximetric Determination

10/22/2019

11. Add about 0.1 g of calcon indicator to 1 of the flasks. Record the initial volume to the nearest 0.01 mL, then titrate to a blue end point. (not to overshoot and compare the colors with reference sample ~~to~~, which has calcon in excess EDTA, to check for the end point).
12. Let the flask stand for 5 minutes after reaching the end point to ensure any  $\text{Ca(OH)}_2$  that precipitates can dissolve and react with EDTA.
13. Continue the titration to a blue, permanent end point (5-10 minutes), if the color of the titrant solution has changed back from blue to purple.
14. Repeat the steps with the remaining samples.
15. Dispose of the EDTA standard solution, the titration solutions & other solution as indicated.



52

## Lab 6

Compleximetric determination

10/22/2019

Partner: Linh Nguyen

Analysis & Results

The analysis and calculations were made in the Excel file included together with this <sup>lab</sup> report. The concentration of each ions in the water sample are as affected in this book.

Conclusion

Two different pH (10 & 12) solutions were prepared to determine the concentration of  $\text{Ca}^{2+}$  &  $\text{Mg}^{2+}$  ions in the water sample. Eriochrome Black T indicator is used to form a red complex in the water sample with either  $\text{Ca}^{2+}$  &  $\text{Mg}^{2+}$  ions. In this pH 10 solutions, the total mol of both ions are determined. At pH 12 meanwhile,  $\text{NaOH}$  was used as indicator, and  $\text{NaOH}$  was added to isolate the  $\text{Mg}^{2+}$  ions and formed  $\text{Mg}(\text{OH})_2$  solid, in order to determine the concentration of  $\text{Ca}^{2+}$  ions only. From here, the concentration of  $\text{Mg}^{2+}$  ions are determined by finding the difference between the average total moles of  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  in the pH 10 solution and each beakers' concentration of  $\text{Ca}^{2+}$ . Then, the concentration of  $\text{MgSO}_4$  and  $\text{CaCl}_2$  in the water sample in mg/L can be determined.