

26

Lab 4 Potentiometric Determination


10/8/2019

Partner: Linh Nguyen

Example Data for Prelab

Fluoride Standards		Cell Potential (mV)	
Concentration (M)	Log[F ⁻] (M)	Average (n = 5)	Std. Dev
0.000001	-6	111	9.8
0.00001	-5	63.8	7.8
0.0001	-4	5.5	6.2
0.001	-3	-52.9	4.4
0.01	-2	-111.1	5.1
0.1	-1	-168.9	3.5

Prelab Graph



$$\%w NaF = \frac{8.5136 \times 10^{-4} g}{0.2215 g} \times 100\% = 0.38\%$$

Sample weighed (g)	0.2215
Volume dissolution (mL)	100
n	5
Average response (mV)	-14.5
Error in response (mV)	5.2
Concentration of F ⁻ (M)	0.000202762
Molar mass of F ⁻ (g/mol)	18.998
Error in molar mass (g/mol)	0.0000005
Mass of F ⁻ (g)	0.000385208
Weight % of F ⁻ (%)	0.17

$$y = -56.646x - 223.69$$

$$56.646x = -223.69 - y \quad \leftarrow \text{avg. response}$$

$$x = \frac{-223.69 - (-14.5)}{56.646}$$

$$x = \log_{10} [F^-]$$

$$\log_{10} [F^-] = \frac{-223.69 + 14.5}{56.646}$$

$$[F^-] = 0.000202762 \text{ M}$$

$$[F^-] = 0.000202762 \text{ M}$$

$$\text{Mass of } F^- = 0.000202762 \text{ mol} \times (0.14) \times$$

$$\left(\frac{18.998 \text{ g}}{\text{mol}} \right) = 0.000385208 \text{ g}$$

$$[NaF] = 0.000202762 \text{ M}$$

$$\text{Mass} = 0.000202762 \text{ mol} \times (0.14) \times (41.988 \text{ g/mol}) \quad \text{Weight \% } F^- = \frac{0.000385208 \text{ g}}{0.2215 \text{ g}} \times 100\% = 0.17\%$$

Potentiometric Determination

10/8/2019

LINEST FUNCTION			
m	-56.6457	-223.6933	b
um	0.8984	3.4987	ub
R ²	0.9990	3.7582	sy
Fisher value	3975.6542	4.0000	df
regression ss	56152.8966	56.4968	residual ss

95% CI	
t table	2.7764
t*um	2.4943
t*ub	9.7140
t* error in average response	14.4375
numerator value	209.1933
error in numerator	17.4012
denominator value	-56.6457
num/denom	-3.6930
error num/denom	-0.3476
Concentration of F- (M)	0.00020
error in concentration (M)	-0.00016

$$x_{\pm ex} = (-223.69 \pm t_{ub}) \div (-14.5 \pm t^*_{\text{error average response}})$$

$$56.646 \pm t_{um}$$

$$e_{num} = \sqrt{(t_{ub})^2 + (t_{ur})^2} = 17.4012 \checkmark$$

$$e_{deno} = 2.4943$$

$$e_x = \sqrt{\left(\frac{17.4012}{(-209.19)}\right)^2 + \left(\frac{2.4943}{-56.6457}\right)^2}$$

$$e_x = \pm 0.3476$$

$$e_y = 0.000202762 (\ln 10)(0.3476)$$

$$e_y = \pm 0.00016$$

$$[F^-] = 0.00020 \pm 0.00016 \text{ M}$$

$$y = 10^x$$

$$\hookrightarrow \frac{e_y}{y} = (\ln 10) e_x$$

$$y = 0.000202762$$

$$e_x = 0.3476$$

28

Lab 4

Partner: Linh Nguyen

Potentiometric Determination

10/8/2019

pH 10/8/19

Actual concentration of stock solution

 $\rightarrow 1015 \text{ ppm} \pm 3 \text{ ppm}$ Two-fold serial dilution

Number of dilutions	1	2	3	4	5
Concentration (ppm)	20.3	10.15	5.075	2.538	1.269

Potentiometric Determination

29

10/8/2019

10/8/2019

Objective

This experiment is conducted to make the students used to using the LabQuest device in order to determine the ionic strength of certain ions. In this lab, we will determine the fluoride potentiometric in different several samples.

Introduction

Some methods learnt from previous lab will be implemented in this experiment, such as ~~10~~-fold ~~10~~ serial dilution. In this lab specifically, a calibration curve is created in order to determine the unknown concentration, based on its potentiometric measurements from the LabQuest device.

Procedure

A) Prepare the standard calibration solutions.

1. All dispensing of unknowns and sample are done by TA.
2. Prepare 20-ppm F^- from the stock solution
(2.00 mL of 1000-ppm stock solution diluted into 100-mL volumetric flask)
3. The dilution is prepared by using the TISAB solution.
4. From that 20-ppm solution, prepare 10-, 5-, 2.5- and 1.25-ppm calibrating fluoride solutions by 2-fold serial dilutions.
5. Use the actual values of ppm (the nominal concentration is just as a guide).

30

Lab 4 Potentiometric Determination
10/8/2019

Partner: Linh Nguyen

First Unknown

Mass of beaker (g)

Mass of beaker + toothpaste (g)

Mass of toothpaste (g)

120.8300 120.3800

120.6624

0.2224

Unknown data

Average cell potential (mV)

11.2

$\log_{10} [F^-] \text{ (ppm)}$

0.536

$[F^-] \text{ (ppm)}$

3.437

weight % of F^- in toothpaste (%)

0.155

error in weight % of F^- in toothpaste, e_w (%)

0.053

Original $[F^-]$ in toothpaste (ppm)

1545

error in original $[F^-]$ in toothpaste, e_e (ppm)

532

Potentiometric Determination

10/8/2019

B) First Unknown: Obtain & Prepare the Toothpaste Sample for Analysis.

- 1) Obtain a toothpaste sample from the instructor.
- 2) A clean, dry and pre-weighed 250-mL beaker should be used and 0.20 to 0.25 g of toothpaste will be added.
- 3) Weigh the beaker with the toothpaste to determine the net weight of the toothpaste sample.

4) Add 50 mL TISAB.

5) Boil the solution gently for 2-3 minutes.

6) Filter the suspension by gravity filtration using #1 or #5 filter paper.

7) After cooling, transfer the solution to a 100-mL volumetric flask and make up to volume with TISAB solution.

C) Second Unknown: Obtain the second F-unknown and note its label.

1) Dispensing of unknowns & stock solutions should be done by the TA only.

2) Use the clean, dry vessel to place the second unknown, which will be analyzed as is (no dilution).

32

Lab 4 Potentiometric Determination 10/8/2019

Partner: Linh Nguyen

pH 0.1

Standard Solution	1	2	3	4	5
Concentration (ppm)					
Potential E(mv), 1		-278.2	-264.8	-249.1	-260.5
Creplicates) 2		-267.6	-264.3	-280.8	-256.0
3			-265.3	-251.9	-256.0
4			-269.9	-251.9	-250.0
5			-263.5	-253.7	-246.1

copied these values
wrongly from
partner.

Unknown #	1st	2nd
Potential E(mv), 1	11.0	
Creplicates) 2	11.4	
3	10.9	
4	11.2	
5	11.6	

no 2nd unknown needed

Standard solution	1	2	3	4	5
Concentration (ppm)	20.3	10.15	5.075	2.538	1.269
Potential E(mv) 1	-33.0	-15.2	-0.15	15.8	39.0
Creplicates) 2	-33.2	-15.1	0.06	15.0	38.9
3	-33.4	-15.9	-0.30	16.3	39.5
4	-33.2	-16.0	-0.18	16.1	40.1
5	-33.3	-15.0	-0.20	16.0	39.7

Lab 4
Potentiometric Determination
10/8/2019

Partner: Linh Nguyen

33

b) Making Potentiometric Measurements with ISE

1) Transfer ~ 20 mL of each standard solution into a clean, dry 50-mL beaker

2) Be sure the solution is well stirred and record ^{to} the cell potential E (measured in mV)

3) The reading will drift for about 1 to 3 mins, wait for the signal to stabilize.

4) Make 4 to 5 replicate measurements for all standards, being careful not to cross-contaminate the solutions.

5) Rinse the electrode with small volumes of the next solution to be tested.

6) Repeat these steps for both unknown sample solutions.

7) Make 4-5 replicate measurements for both unknown sample solutions.

8) Make sure the voltage of unknown samples falls within the voltage range of the standards; or else, preparation & measurement of additional calibrating solutions may be needed, as extrapolation is not a good idea.

9) Dispose all solutions as indicated by the instructor. Rinse all glassware and clean the work area. Return all glassware.

34


Lab 4 Potentiometric Determination
10/8/2019

Partner: Linh Nguyen
10/8/2019

Analysis & Results

The analysis and extra calculations were done in excel file and it was included together with this lab's summary report. The standard calibration curve is as shown below. It follows the normal Nernstian operation, with a slope within -59 ± 4 mV. The blue point is the fluoride concentration of toothpaste sample's result.

Calibration curve of Fluoride Ions



Conclusion

The toothpaste sample's cell potential measurement falls within the standard calibration curve, in which the concentration of fluoride ions in the toothpaste can be determined, which was 1545 ± 532 ppm.