



Exp. No. 5	Experiment/Subject CHEM 438 (POF)	Date 10/11/2021
Name Abdul Fayed	Lab Partner Justin M	Locker/Desk No. Course & Section No. 021L


PRE-LAB Questions

- $\lambda_{ex} = 340\text{nm}$ (peak @ 3400\AA) \rightarrow no peak for toluene @ this λ .
 $\lambda_{em} = 460\text{nm}$ (peak @ 4600\AA)
- $\epsilon = 9750\text{M}^{-1}\text{cm}^{-1}$ (from the graph @ $\lambda_{obs} = 340\text{nm}$)
 $\Omega = 2$, $P_0 = 0.5$, $b = 1\text{cm}$



(The specified conc. range) (10-fold dilution range)
 \Rightarrow should dilute the quinine sample to get linear range of calib. curve.

- Tonic water should be diluted. 40-dilution factor ~~should~~ ^{can} be used.



- Take average of blank values then use the standard deviation to calculate the limit of detection.

Exp. No. <u>65</u>	Experiment/Subject <u>CHEM 438 (D&F)</u>	Date <u>10/11/2021</u>
Name <u>Abdul Fayeed</u>	Lab Partner <u>Justin M</u>	Locker/Desk No. <u>0211</u>

Objective

Purpose of this lab is to determine the concentration of quinine in tonic water by fluorescence.

Introduction

Fluorescence is ~~effective~~ selective to one or 2 compounds in a mixture. Tonic water has both quinine and benzoate. Appropriate λ_{ex} and λ_{em} have to be picked to ensure that only one compound can be fluoresce at a time. Fluorescence spectrum also has its limit where it reaches a maximum fluorescence then drop down as concentration increases.

Procedure

1. Determine λ_{ex} and λ_{em} .
 $\lambda_{ex} = 340\text{nm}$, $\lambda_{em} = 460\text{nm}$
2. Set non-scanned (excitation) wavelength to appropriate λ_{ex} value.
3. Measure fluorescence excitation spectrum of quinine once, using the highest concentration, using λ_{em} to scan the excitation wavelength.
4. Measure the fluorescence emission spectra of quinine by fixing λ_{ex} and scanning over emission wavelength.
5. Prepare a set of fluorescence calibration standards to determine the diluted tonic

mass of quinine sulfate (g) =
stock concentration =

sample	Vol of stock (mL)	conc (M)	F_1
1			
2			
3			
4			
5			
Tonic			
Blank	0	0	


sample	F_2	F_3
1		
2		
3		
4		
5		
Tonic		
Blank		

Quinine 'unknown'. Dilute solution w/ 0.05M H₂SO₄.

6. Make measurement to determine the LOD of fluorescence. (10 times)

Conclusions

Quinine can be detected and measured based on its concentration by fluorescence. Fluorescence is selective to only one or 2 compounds, so it can save time to separate the compound from the mixture.


Signature 	Date <u>10/11/2021</u>	Witness/TA	Date
---	------------------------	------------	------

Exp. No. 5	Experiment/Subject CHEM 438 (DGF)	Date 10/11/2021
Name Abdul Fayed	Lab Partner Justin M	Locker/ Desk No.
		Course & Section No. 0211

LOD data

Run #

Response

Signature 	Date 10/11/2021	Witness/TA	Date
---	-----------------	------------	------