


Sif


Lab 7
Hypochlorite Ion in Bleach
10/29/2019

Partner: Linh Nguyen
10/29/2019

Relab Assignment




$$0.1358\text{g} \left(\frac{1\text{mol}}{214.00\text{g}} \right) = 6.34579 \times 10^{-4} = n_{\text{IO}_3^-} = n_{\text{IO}_3^-}$$




$$n_{\text{I}_3^-} = 1.903737 \times 10^{-3} \text{ mol} \quad (a)$$

$$2\text{g} \left(\frac{1\text{mol}}{166.0028\text{g}} \right) = 0.01205 \text{ mol KI} = 0.01205 \text{ mol I}^-$$




$$n_{\text{I}_3^-} = 4.51875 \times 10^{-3} \text{ mol} \quad (b)$$

$$0.002\text{L}(6\text{M}) = 0.012 \text{ mol} = n_{\text{HCl}} = n_{\text{H}^+}$$



$$n_{\text{I}_3^-} = 6 \times 10^{-3} \text{ mol} \quad (c)$$

Since $n_{\text{I}_3^-} (a) < (b) < (c)$, IO_3^- is the limiting reagent.



$$n_{\text{S}_2\text{O}_3^{2-}} = 2(1.903737 \times 10^{-3} \text{ mol}) = 3.807474 \times 10^{-3}$$

$$[\text{S}_2\text{O}_3^{2-}] = \frac{3.807474 \times 10^{-3} \text{ mol}}{0.03194\text{L}} = 0.12 \text{ M}$$

$$\frac{M_1 V_1}{M_2 V_2} = 2$$

$$0.12\text{M}(21.33\text{mL}) = 2 \times M_2(25.00\text{mL})$$

$$M_2(25.00\text{mL})$$

$$M_2 = 0.05 \text{ M} = [\text{I}_3^-]$$



10/29/2019

Objective

The purpose of this experiment is to determine the concentration of hypochlorite ions in bleach by titration of thiosulfate ions solution into the bleach solution.

Introduction

The experiment is divided into two parts, where the first part involves the determination of concentration of thiosulfate solution, to be used in the second part, to then determine the concentration of hypochlorite ions in bleach, via titrations. Starch is used as an indicator, which will change from blue to black colored complex with iodine. For iodimetry (where standard iodine solution is used as titrant for oxidizable titrant), the end point corresponds to a sudden color change to blue. Meanwhile for iodometry (where excess iodine is used to reduce a chemical species while simultaneously forming iodine), the end point corresponds to a sudden loss of blue color in the solution.

Procedure

A) Preparation of Standard Sodium Thiosulfate solution

1. Prepare a 0.1M thiosulfate solution (dissolve ~12g of sodium thiosulfate pentahydrate & ~0.05g of sodium carbonate in ~250 mL of water). Mix very well to dissolve. Transfer the solution to a clean 500 mL volumetric flask. Dilute to the mark & mix well. (serve as stock thiosulfate solution).
2. Fill a 50 mL burette with thiosulfate solution. Record the initial volume.

56

Lab. paper 10/29/2019
Hypochlorite ion in Bleach

Partner: Linh N. Nguyen

10/29/2019

Part A

Erlenmeyer flask #	1	2	3
Mass of KI (g)	1.8801	2.1689	2.1266
Beakers #	1	2	3
Mass of KIO_3 solid (g)	0.1421	0.1249	0.1376
Initial buret reading (mL)	0.29	0.23	0.82
Final buret reading (mL)	42.35	36.65	42.18
Titrated volume (mL)	42.06	36.42	41.36

Sodium Thiosulfate standardization

Erlenmeyer flask #	1	2	3
Mass of KI solid (g)	1.8801	2.1689	2.1266
Mass of KIO_3 solid (g)	0.1421	0.1249	0.1376
Initial buret reading (mL)	0.29	0.23	0.82
Final buret reading (mL)	42.35	36.65	42.18
Titrated volume (mL)	42.06	36.42	41.36
Molecular weight of KIO_3 (g/mol)	214.001		
Mole of IO_3^- ions (mol)	6.640E-04	5.836E-04	6.430E-04
Mole of I_3^- ions (mol)	1.992E-03	1.751E-03	1.929E-03
Mol of $\text{S}_2\text{O}_3^{2-}$ ions (mol)	3.984E-03	3.502E-03	3.858E-03
Molarity of $\text{Na}_2\text{S}_2\text{O}_3$ solution (M)	0.09472	0.09615	0.09328
Average molarity (M)	0.09472		
Standard deviation in molarity (M)	0.00144		
RSD (%)	1.5		



Hypochlorite Ion in Bleach

10/29/2019

3. Weigh $\sim 2g$ of reagent-grade KI solid into each of 3 numbered 125 mL Erlenmeyer flasks. Add ~ 25 mL of water & ~ 2 mL of 6 M HCl to each flask. Swirl the solution to dissolve the KI.
4. Number 3 small beakers (with the same numbers used for Erlenmeyer flasks).
5. Into each beaker, weigh between 0.12 & 0.15 g of KIO_3 . Record the mass of KIO_3 in each numbered container/beaker.
6. Add KIO_3 in one of the numbered containers to the solution in the identically numbered Erlenmeyer flask. Swirl the solution to dissolve the potassium KIO_3 . Brown color of I_3^- should be apparent in the solution.
7. After recording the initial volume, titrate the solution in the Erlenmeyer flask with the thiosulfate solution in the burette until the I_3^- color has become noticeably less intense. Add ~ 2 mL of starch indicator & continue the titration until the dark-blue color of starch-triiodide complex disappears. Record end-point volume to the nearest 0.01 mL.
8. Individually repeat steps 6 & 7 with the remaining 2 KIO_3 samples. Generate iodine/triiodide just before titrating with thiosulfate to minimize loss of triiodide reactant & error from the loss.

58

Lab 7

Hypochlorite Ion in Bleach

10/29/2019

Partner: Linh Nguyen

Part B

Erlenmeyer flasks #	1	2	3
Mass of KI (g)	1.5865	1.9429	1.7282
Initial buret reading (mL)	42.18	14.93	18.95
Final buret reading (mL)	45.90	18.95	23.49
Titrated volume (mL)	3.72	4.02	4.54

Determination of Hypochlorite in Bleach

Erlenmeyer flask	1	2	3
Mass of KI solid (g)	1.5865	1.9429	1.7282
Initial buret reading (mL)	42.18	14.93	18.95
Final buret reading (mL)	45.90	18.95	23.49
Titrated volume (mL)	3.72	4.02	4.54
Average molarity $\text{Na}_2\text{S}_2\text{O}_3$ solution (M)	0.09472		
Volume of NaOCl (mL)	2.00		
Mole of $\text{S}_2\text{O}_3^{2-}$ ions (mol)	3.523E-04	3.808E-04	4.300E-04
Mole of I_3^- ions (mol)	1.762E-04	1.904E-04	2.150E-04
Mole of OCI^- (mol)	1.762E-04	1.904E-04	2.150E-04
Molarity of NaOCl (M)	0.08809	0.09519	0.1075
Average molarity NaOCl (M)	0.09693		
Standard deviation molarity NaOCl (M)	0.00982		
RSD (%)	10.1		



10/29/2019

B) Iodometric Determination of Hypochlorite

1. Add ~ 50 mL of DI water & between 1.5 & 2.0g of KI to each of 3 250 mL Erlenmeyer flasks. Stir the solutions until the ~~pot~~ KI in each flask is dissolved.
2. Add ~ 10 mL of 3M H_2SO_4 to each flask.
3. Use a pipet to deliver 2.00 mL of liquid bleach into one of the flasks, & stir by swirling the solution. Fill a 50 mL burette with standard sodium thiosulfate solution, record the initial volume, & use the thiosulfate solution to titrate the I_3^- formed during the reaction of bleach with iodide. When the I_3^- color starts to fade, add 2 mL of starch solution & continue the titration to the end point. Record the end-point volume to the nearest 0.01 mL.
4. Individually repeat step 3 with the solutions in each of the other 2 Erlenmeyer flasks. Generate iodine just before the titration with thiosulfate.
5. Dispose of the titration solutions, thiosulfate standard & other solutions as indicated. Clean & dry the lab work area.

60

Lab 7

Hypochlorite Ion in Bleach

10/29/2019

Partner: Linh Nguyen

Analysis & Results

The calculations for the determination of molarity of thiosulfate solution and the liquid bleach are done in Excel, as tabulated in the table shown attached in this book. The formula used were as stated in the summary report.

Conclusion

In order to determine the concentration of liquid bleach, the thiosulfate solution has to be prepared. The average molarity of thiosulfate solution is found to be $0.09422M$ with a standard deviation of $0.00144M$. By using this thiosulfate solution, it is titrated into the liquid bleach solution and the molarity of the bleach is found to be $0.09693M$ with a standard deviation of $0.00982M$.